



Kids' Medical Care
Diana McLaughlin, MD
2336 Immokalee Road
Naples, FL 34110
239-591-8481
239-596-0212 (FAX)

Name: _____
MR#: _____
DOB: _____
Date: _____

NEBULIZER MACHINE INSTRUCTIONS

You may watch a video on how to use a nebulizer machine by going to www.pediatricsupplier.com.

1. Wash your hands before handling the nebulizer unit and tubing.
2. Plug the unit into a properly working electrical outlet.
3. Place the unit on a firm surface such as a table top.
4. Do not cover vents on the unit.
5. Connect one end of the tubing to the unit and the other end to the bottom of the medication cup.
6. Add medication to the cup along the side not down the center.
7. Replace the top of the medication cup and attach either a mask or mouthpiece. Place the mask or mouthpiece on the child.
8. Turn the power switch on.
9. Hold the medication cup upright. Tilting or holding cup sideways will cause the meds to spill out or not mist correctly.
10. Breathe comfortably, slowly and deeply. Breathing rapidly may cause hyperventilation. It is acceptable if the child is crying. The child may receive a nebulizer treatment while sleeping.
11. The nebulizer treatment is finished when you hear a sputtering sound or there is no mist. It is acceptable if there are a few drops of meds still left. The treatment may take from 5 to 15 minutes.
12. Turn the unit off and unplug the unit from the wall.
13. Wash mask or mouthpiece and disassembled medication cup with mild soap and warm water. Do not wash tubing.
14. Dry with lint free cloth or allow to air dry. Do not expose to direct sunlight.
15. Change out the filter once a week during regular use.
16. You may call the office if you have any questions. If the unit malfunctions, please call Pediatric Suppliers at 1-866-380-3484.

I, certify that I have received and understand the above instructions.

Print Name

Signature

Date

INSTRUCTOR